

## GENERAL INFORMATION

**Product name :** Jointed adapter (male) with PTFE key stopcock**Description :** Made of borosilicate glass LBG 3.3. 9 mm hose connector for 8 mm Øint tubing. Angled or straight

## TECHNICAL DATA

reference	joint ST/NS	description	pcs/pack
ADFL-001-001	14/23	angled connection	1
ADFL-002-001	19/26	angled connection	1
ADFL-003-001	24/29	angled connection	1
ADFL-004-001	29/32	angled connection	1
ADFL-005-001	14/23	straight connection	1
ADFL-006-001	19/26	straight connection	1
ADFL-007-001	24/29	straight connection	1
ADFL-008-001	29/32	straight connection	1

## PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
ADFL-001-001	0,4	0,15	70172000	08434868088214
ADFL-002-001	0,4	0,15	70172000	08434868088221
ADFL-003-001	0,4	0,15	70172000	08434868088238
ADFL-004-001	0,4	0,15	70172000	08434868088245
ADFL-005-001	0,4	0,15	70172000	08434868088252
ADFL-006-001	0,4	0,15	70172000	08434868088269
ADFL-007-001	0,4	0,15	70172000	08434868088276
ADFL-008-001	0,4	0,15	70172000	08434868088283

PRODUCT PHOTO



## INFORMACIÓN GENERAL

**Nombre del producto : Adaptador esmerilado macho con llave de PTFE****Descripción :** De vidrio borosilicato LBG 3.3. Con oliva de diámetro 9 mm para tubos de Øint 8 mm. Curvado o recto

## DATOS TÉCNICOS

referencia	esmerilado	descripción	unidades por ref.
ADFL-001-001	14/23	curvado	1
ADFL-002-001	19/26	curvado	1
ADFL-003-001	24/29	curvado	1
ADFL-004-001	29/32	curvado	1
ADFL-005-001	14/23	recto	1
ADFL-006-001	19/26	recto	1
ADFL-007-001	24/29	recto	1
ADFL-008-001	29/32	recto	1

## EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
ADFL-001-001	0,4	0,15	70172000	08434868088214
ADFL-002-001	0,4	0,15	70172000	08434868088221
ADFL-003-001	0,4	0,15	70172000	08434868088238
ADFL-004-001	0,4	0,15	70172000	08434868088245
ADFL-005-001	0,4	0,15	70172000	08434868088252
ADFL-006-001	0,4	0,15	70172000	08434868088269
ADFL-007-001	0,4	0,15	70172000	08434868088276
ADFL-008-001	0,4	0,15	70172000	08434868088283

FOTO DEL PRODUCTO



## INFORMATIONS GÉNÉRALES

**Nom produit : Raccord rodé, mâle, avec clé en PTFE****Description :** En verre borosilicaté LBG 3.3. Avec olive de diamètre 9 mm pour tubes de Øint 8 mm. Coudé ou droit

## DONNÉES TECHNIQUES

référence	rodage	description	unités par ref.
ADFL-001-001	14/23	coudé	1
ADFL-002-001	19/26	coudé	1
ADFL-003-001	24/29	coudé	1
ADFL-004-001	29/32	coudé	1
ADFL-005-001	14/23	droit	1
ADFL-006-001	19/26	droit	1
ADFL-007-001	24/29	droit	1
ADFL-008-001	29/32	droit	1

## EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
ADFL-001-001	0,4	0,15	70172000	08434868088214
ADFL-002-001	0,4	0,15	70172000	08434868088221
ADFL-003-001	0,4	0,15	70172000	08434868088238
ADFL-004-001	0,4	0,15	70172000	08434868088245
ADFL-005-001	0,4	0,15	70172000	08434868088252
ADFL-006-001	0,4	0,15	70172000	08434868088269
ADFL-007-001	0,4	0,15	70172000	08434868088276
ADFL-008-001	0,4	0,15	70172000	08434868088283



## INFORMAZIONE GENERALE

**Nome del prodotto :** Raccordo smerigliato con rubinetto in PTFE**Descrizione :** In vetro borosilicato LBG 3.3. Il raccordo è per tubi con Ø<sub>int</sub> da 8 mm. Curvo o dritto

## DATI TECNICI

referenza	smerigliato	descrizione	unità per ref.
ADFL-001-001	14/23	curvo	1
ADFL-002-001	19/26	curvo	1
ADFL-003-001	24/29	curvo	1
ADFL-004-001	29/32	curvo	1
ADFL-005-001	14/23	dritto	1
ADFL-006-001	19/26	dritto	1
ADFL-007-001	24/29	dritto	1
ADFL-008-001	29/32	dritto	1

## IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
ADFL-001-001	0,4	0,15	70172000	08434868088214
ADFL-002-001	0,4	0,15	70172000	08434868088221
ADFL-003-001	0,4	0,15	70172000	08434868088238
ADFL-004-001	0,4	0,15	70172000	08434868088245
ADFL-005-001	0,4	0,15	70172000	08434868088252
ADFL-006-001	0,4	0,15	70172000	08434868088269
ADFL-007-001	0,4	0,15	70172000	08434868088276
ADFL-008-001	0,4	0,15	70172000	08434868088283

