

GENERAL INFORMATION

Product name : Latex tubing**Description :** Amber colour. For lab applications where liquids must be transferred

TECHNICAL DATA

reference	Øinside (mm)	Øouter (mm)	type	pcs/pack
LTXT-004-005	4.5	9	for 4/5 hose	5
LTXT-007-005	7	10	for 8/9 hose	5
LTXT-008-005	8	12	for 10 hose	5
LTXT-010-005	10	14	for 12/13 hose	5
LTXT-012-005	12	18	for 15/16 hose	5
LTXT-014-005	14	20	for 16 hose	5
LTXT-004-025	4.5	9	for 4/5 hose	25
LTXT-007-025	7	10	for 8/9 hose	25
LTXT-008-025	8	12	for 10 hose	25
LTXT-010-030	10	14	for 12/13 hose	30
LTXT-012-025	12	18	for 15/16 hose	25
LTXT-014-020	14	20	for 16 hose	20

PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
LTXT-004-005	0,01	0,001	40029100	08434868028203
LTXT-007-005	0,75	0,25	40029100	08434868028227
LTXT-008-005	0,01	0,001	40029100	08434868028241
LTXT-010-005	0,01	0,001	40029100	08434868028265
LTXT-012-005	0,01	0,001	40029100	08434868028289
LTXT-014-005	0,01	0,001	40029100	08434868028302
LTXT-004-025	4	1,5	40029100	08434868028210
LTXT-007-025	3	1,1	40029100	08434868028234
LTXT-008-025	0,01	0,001	40029100	08434868028258
LTXT-010-030	0,01	0,001	40029100	08434868028272
LTXT-012-025	0,01	0,001	40029100	08434868028296
LTXT-014-020	0,01	0,001	40029100	08434868028319

PRODUCT PHOTO



INFORMACIÓN GENERAL

Nombre del producto : Tubería de látex**Descripción :** Color ámbar. Para refrigerantes y otros montajes en los que deba hacerse circular líquidos

DATOS TÉCNICOS

referencia	Øinterior (mm)	Øexterior(mm)	tipo	unidades por ref.
LTXT-004-005	4.5	9	para olivas 4/5	5
LTXT-007-005	7	10	para olivas 8/9	5
LTXT-008-005	8	12	para olivas 10	5
LTXT-010-005	10	14	para olivas 12/13	5
LTXT-012-005	12	18	para olivas 15/16	5
LTXT-014-005	14	20	para olivas 16	5
LTXT-004-025	4.5	9	para olivas 4/5	25
LTXT-007-025	7	10	para olivas 8/9	25
LTXT-008-025	8	12	para olivas 10	25
LTXT-010-030	10	14	para olivas 12/13	30
LTXT-012-025	12	18	para olivas 15/16	25
LTXT-014-020	14	20	para olivas 16	20

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
LTXT-004-005	0,01	0,001	40029100	08434868028203
LTXT-007-005	0,75	0,25	40029100	08434868028227
LTXT-008-005	0,01	0,001	40029100	08434868028241
LTXT-010-005	0,01	0,001	40029100	08434868028265
LTXT-012-005	0,01	0,001	40029100	08434868028289
LTXT-014-005	0,01	0,001	40029100	08434868028302
LTXT-004-025	4	1,5	40029100	08434868028210
LTXT-007-025	3	1,1	40029100	08434868028234
LTXT-008-025	0,01	0,001	40029100	08434868028258
LTXT-010-030	0,01	0,001	40029100	08434868028272
LTXT-012-025	0,01	0,001	40029100	08434868028296
LTXT-014-020	0,01	0,001	40029100	08434868028319



INFORMATIONS GÉNÉRALES

Nom produit : Tuyau en latex**Description :** Couleur jaune. Pour réfrigérant et autre montage nécessitant la circulation de liquide

DONNÉES TECHNIQUES

référence	Øintérieur (mm)	Øextérieur (mm)	type	unités par ref.
LTXT-004-005	4.5	9	pour olive 4/5	5
LTXT-007-005	7	10	pour olive 8/9	5
LTXT-008-005	8	12	pour olive 10	5
LTXT-010-005	10	14	pour olive 12/13	5
LTXT-012-005	12	18	pour olive 15/16	5
LTXT-014-005	14	20	pour olive 16	5
LTXT-004-025	4.5	9	pour olive 4/5	25
LTXT-007-025	7	10	pour olive 8/9	25
LTXT-008-025	8	12	pour olive 10	25
LTXT-010-030	10	14	pour olive 12/13	30
LTXT-012-025	12	18	pour olive 15/16	25
LTXT-014-020	14	20	pour olive 16	20

EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
LTXT-004-005	0,01	0,001	40029100	08434868028203
LTXT-007-005	0,75	0,25	40029100	08434868028227
LTXT-008-005	0,01	0,001	40029100	08434868028241
LTXT-010-005	0,01	0,001	40029100	08434868028265
LTXT-012-005	0,01	0,001	40029100	08434868028289
LTXT-014-005	0,01	0,001	40029100	08434868028302
LTXT-004-025	4	1,5	40029100	08434868028210
LTXT-007-025	3	1,1	40029100	08434868028234
LTXT-008-025	0,01	0,001	40029100	08434868028258
LTXT-010-030	0,01	0,001	40029100	08434868028272
LTXT-012-025	0,01	0,001	40029100	08434868028296
LTXT-014-020	0,01	0,001	40029100	08434868028319



INFORMAZIONE GENERALE

Nome del prodotto : Tubo in lattice**Descrizione :** Colore ambrato. Per refrigeranti e altri sistemi in cui devono circolare liquidi

DATI TECNICI

referenza	Øinterno (mm)	Øesterno (mm)	tipo	unità per ref.
LTXT-004-005	4.5	9	per portagomme 4/5	5
LTXT-007-005	7	10	per portagomme 8/9	5
LTXT-008-005	8	12	per portagomme 10	5
LTXT-010-005	10	14	per portagomme 12/13	5
LTXT-012-005	12	18	per portagomme 15/16	5
LTXT-014-005	14	20	per portagomme 16	5
LTXT-004-025	4.5	9	per portagomme 4/5	25
LTXT-007-025	7	10	per portagomme 8/9	25
LTXT-008-025	8	12	per portagomme 10	25
LTXT-010-030	10	14	per portagomme 12/13	30
LTXT-012-025	12	18	per portagomme 15/16	25
LTXT-014-020	14	20	per portagomme 16	20

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
LTXT-004-005	0,01	0,001	40029100	08434868028203
LTXT-007-005	0,75	0,25	40029100	08434868028227
LTXT-008-005	0,01	0,001	40029100	08434868028241
LTXT-010-005	0,01	0,001	40029100	08434868028265
LTXT-012-005	0,01	0,001	40029100	08434868028289
LTXT-014-005	0,01	0,001	40029100	08434868028302
LTXT-004-025	4	1,5	40029100	08434868028210
LTXT-007-025	3	1,1	40029100	08434868028234
LTXT-008-025	0,01	0,001	40029100	08434868028258
LTXT-010-030	0,01	0,001	40029100	08434868028272
LTXT-012-025	0,01	0,001	40029100	08434868028296
LTXT-014-020	0,01	0,001	40029100	08434868028319

