


qPCR 8-TUBE STRIPS
PCRL-384-100
 astik'S
 Date: 08-03-2023

Material: Polypropylene USP VI
 Colour: Natural, transparent
 Capacity: 0,1 ml
 Features: Alphanumeric code. Full skirted. Compatible with ABI & Bio-Rad.
 Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 8-TUBE STRIPS
PCRL-W84-100

astik's

Date: 08-03-2023

Material: Polypropylene USP VI

Colour: white, opaque

Capacity: <0,1 ml

Features: Alphanumeric code. Full skirted. Compatible with Roche & Bio-Rad.

Human DNA free and certified DNase, RNase, and PCR inhibitors free