

	Capacity	Standard Line	Premium Line	ULR Premium Line	Filtered Premium Line
<b>EASY 80+</b>					
EASC-10P-001	100 µl	TIPP-200-***	TICD-200-***	TICU-200-***	TICF-100-***, TICF-200-***
EASC-20P-001	200 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASC-50P-001	500 µl	TIPP-1K0-***, TIPP-1K1-***	TICD-1K0-***	TICU-1K0-***	TICF-1K0-***
EASC-1KP-001	1000 µl	TIPP-1K0-***, TIPP-1K1-***	TICD-1K0-***	TICU-1K0-***	TICF-1K0-***
<b>EASY 40+</b>					
EASD-P02-001	0,1 - 2,5 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-P10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-P20-001	2 - 20 µl	-	-	-	-
EASD-P50-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-100-***
EASD-10P-001	10 - 100 µl	TIPP-200-***	TICD-200-***	TICU-200-***	TICF-100-***, TICF-200-***
EASD-20P-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-1KP-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TICD-1K0-***	TICU-1K0-***	TICF-1K0-***
EASD-5KP-001	500 - 5000 µl	-	TIAD-5K0-***	-	-
EASD-10K-001	1000 - 10000 µl	-	TIAD-10K-***	-	-
<b>EASY 40 ELITE</b>					
EASD-E02-001	0,1 - 2,5 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-E10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-E20-001	2 - 20 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-020-***
EASD-E50-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-100-***
EASD-10E-001	10 - 100 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-20E-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-1KE-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TIAD-1K0-***	TICU-200-***	TICF-1K0-***
EASD-5KE-001	1000 - 5000 µl	-	TIAD-5K0-***	-	-
<b>EASY 40 Xpert</b>					
EASD-X02-001	0,1 - 2,5 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-X10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-X20-001	2 - 20 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-020-***
EASD-X50-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-100-***
EASD-10X-001	10 - 100 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-20X-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-30X-001	30 - 300 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-1KX-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TIAD-1K0-***	TICU-1K0-***	TICF-1K0-***
EASD-5KX-001	1000 - 5000 µl	-	TIAD-5K0-***, TIAD-5K1-***	-	-
EASD-1KX-001	1000 - 10000 µl	-	TIAD-18K-100	-	-
<b>EASY 40 Xpert 8-channels</b>					
EASD-G10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-G20-001	2 - 20 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-020-***
EASD-G50-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-100-***
EASD-10G-001	10 - 100 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-20G-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
<b>EASY 40 Xpert 12-channels</b>					
EASD-T10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASD-T20-001	2 - 20 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-020-***
EASD-T50-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-100-***
EASD-10T-001	10 - 100 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-20T-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASD-30T-001	30 - 300 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
<b>EASY W</b>					
EASW-010-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASW-050-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-100-***
EASW-300-001	30 - 300 µl	TIPP-200-***, TIPP-201-***	TIAD-300-***, TICD-200-***	TICU-200-***	TICF-200-***
EASW-1K0-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TICD-1K0-***	TICU-1K0-***	TICF-1K0-***
<b>DLAB HiPette</b>					
EASH-H02-001	0,1 - 2,5 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASH-H10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TICD-010-***	TICU-010-***	TICF-010-***
EASH-H20-001	2 - 20 µl	-	-	-	-
EASH-H10-001	10 - 100 µl	TIPP-200-***, TIPP-201-***	TICD-200-***	TICU-200-***	TICF-100-***
EASH-H20-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASH-H30-001	30 - 300 µl	TIPP-200-***, TIPP-201-***	TICD-200-***, TIAD-300-***	TICU-200-***	TICF-200-***
EASH-1KH-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TICD-1K0-***	TICU-1K0-***	TICF-1K0-***
EASH-5KH-001	1000 - 5000 µl	-	TIAD-5K0-***	-	-
EASH-10K-001	1000 - 10000 µl	-	TIAD-10K-***	-	-